## CLASS-10 PHYSICAL SCIENCE PERIOD PLANS

## **CHAPTER: 01 – HEAT**

## **PERIOD PLAN-05**: Calorie meter – structure – working

Content Analysis	Class Room Environment	Teaching Learning Material
Calorimeter: A device used to measure the quantity of heat.  Description: Calorimeter consists of a thin cylindrical copper vessel. The copper vessel is placed in a wooden box. The gap between vessel and wooden box is filled with insulating material like wool/fur. The wooden box contains a wooden lid. The lid has a slot that a stirrer is to be immersed through it into the copper vessel. Also the lid has a hole to immerse laboratory thermometer. (It may also have an holder for thermometer at a side of the wooden box.)	Conversation: About the description of calorimeter by showing the parts using calorimeter.  STIRRER  STIRRER  THERMOMETER  THERMOS FLASK  BEAKER  SOLUTION  CORK/WOODEN BLOCK	Calorimeter model and photos
Working procedure: First we have to find the mass of the calorimeter $(m_1)$ . Fill nearly half of the calorimeter with water and find the mass of calorimeter with water ( $m_2$ ). Measure the initial temperature with laboratory thermometer $(T_1^{\circ}C)$ . This is the temperature of both water and also calorimeter.  Take a few lead shots and place them in hot water. Heat them up to a temperature $100^{\circ}C$ . So measure the temperature of lead shots $(T_2^{\circ}C)$ . Transfer the lead shots into calorimeter quickly with minimum loss of heat. Stir the mixture well. Note the final temperature $(T_3^{\circ}C)$ . Measure the final mass of calorimeter along with water and lead shots $(m_3)$ . Heat $(Q) = m.s.\Delta T$ According to the method of mixtures: Heat lost by the solid = Heat gained by calorimeter + Heat gained by water $(m_3-m_2).S_1.(T_2-T_3) = m_1.S_c.(T_3-T_1) + (m_2-m_1).S_w.(T_3-T_1) $ $S_l = \frac{[m_1 S_c + (m_2-m_1)S_w][T_3-T_1]}{(m_3-m_2)(T_2-T_3)}$ This way we can find the specific heat of a solid. $S_w = 1 \text{ cal/gm -}^{\circ}C$ $S_c = 0.095 \text{ cal/gm -}^{\circ}C$	Explanation: Explains the procedure for determination of specific heat of a solid using calorimeter.	AV-Clip

NAGA MURTHY- 9441786635

Contact at: nagamurthysir@gmail.com
Visit at: ignitephysics.weebly.com